

What conditions can be treated?

- Acne
- Enlarged Pores
- Rosacea
- Redness
- Purple Pigment
- Facial Contour
- Uneven Pigment
- Sun Damage
- Wrinkles
- Fine Lines
- Scars
- Sagging Skin

Combination Therapy?


For accelerated ageing or long standing/difficult conditions your therapist may recommend that you combine a course of Soli-Tone 2500 treatments with either surface skin peels, VPL or micro-dermabrasion to ensure your receive rapid and intensive results. If you do require combination therapy your therapist will assist you in planning and scheduling your individual course of treatments.

How soon can I Start!

The most difficult thing is getting an available appointment! Firstly, however, you need to start with a quick consultation to assess your skin and your needs. Consultation is free and takes about 10-20 minutes. Then prepare yourself for the biggest and best changes you will ever see in your skin.

The Rules?

We are skin fanatics. We are also fanatical about results. So in combination with your treatments you **MUST** use prescribed skin care that will compliment your treatments, speed up the results you want, and most importantly maintain your new younger looking skin. Your therapist will recommend suitable products and provide you with samples to try.



Soli-Tone 2500 uses 4 specific wavelengths for instant access to

Nothing else even comes close.
Soli-Tone 2500

Exclusive to: Absolute Beauty

Ph 6331 0449



What would you pay to look years & years & years younger?

well with the Soli-Tone 2500 - NOT MUCH!

Face lift, thousands. Botox hundreds. Now EVERYONE can afford to look young, all the time and without the pain involved with surgery or Botox.

NO Downtime and RESULTS in weeks.

| | |
|--|------------------|
| LumiFacial® | \$140.00 |
| For Acne, Blemishes, Uneven Pigmentation, Redness, Rosacea, Photo-aging | |
| Package 8 LumiFacials® | \$980.00 |
| LumiFacial® & LumiLift® | \$190.00 |
| Includes LumiFacial® for above conditions and LumiLift® for <u>lifting and tightening</u> sagging muscles of the eyes, face and décolletage. | |
| Package 8 LumiFacial® & LumiLift® | \$1330.00 |
| LumiFacial® Add On Treatment | \$95.00 |
| <u>Add</u> to your favourite Absolute Beauty Facial, Peel or Micro-dermabrasion treatment for accelerated results. | |
| Monthly Maintenance (4 weekly) LumiFacial® | \$110.00 |
| <u>Only</u> available to clients who have <i>previously purchased a package</i> of LumiFacials®. | |
| Monthly Maintenance (4 weekly) LumiFacial® & LumiLift® | \$140.00 |
| <u>Only</u> available to clients who have <i>previously purchased a package</i> of LumiFacial® & LumiLift® Treatments. | |

Soli-Tone 2500 exclusive to:



Simply the best in beauty.

Absolute Beauty

175 Elphin Road, Newstead.
(Next to Elphin Continental Cakes)

www.absolutebeauty.com.au

Ph **6331 0449**



Soli-Tone 2500


Photo Biostimulation system.

The most advanced non-surgical face lifting treatment available TODAY.

Take years off your face in weeks. NOW you can wind back the clock.




Before



After

Acute Acne

Before and after 4 weeks of combination therapy. (4 x Peels, 9 x LumiFacial) 13 Treatments all up. Now free from acne, 12 months on. Previously tried medication and other therapy with no results.



Before



After

Mother & Daughter? NO.

Before and after 8 weeks of combination therapy. (6 x LumiFacial/Lift with Pre-Micro Treatment) Looks at least 10 Years Younger

Due to varying print quality these pictures may seem less amazing than the real results are! Please ask to see high resolution graphics in salon. You will be amazed by the results we can achieve for you.

Photobiostimulation system - a uniquely advanced technology for skin rejuvenation. The biggest advance in non-surgical face lifting.

Take YEARS off your face in WEEKS.

Using technology produced by NASA and unlike anything else available on the market today, we can now REVERSE (yes that is reverse) the ageing of your skin and take years off your face in a matter of weeks. The Soli-Tone 2500 is a unique skin rejuvenation device that uses light therapy and special currents to turn back the clock on ageing, sun-damaged skin, blemishes, redness, and uneven pigmentation. This multi-platform system features two treatments: the Lumilift® and the LumiFacial®. Simultaneous application of these technologies mutually enhances their individual effectiveness, promoting cellular regeneration from within. **Results are immediate and visible.**



What is the treatment and how does the treatment work?

The Soli-Tone 2500 is registered with the FDA and uses three technologies: pure spectrum Photo-Energy (LED), Impulse micro currents and ultrasound. LED (Light Emitting Diode) Therapy was initially developed by NASA and has been used to treat serious burns, crush injuries, diabetic ulcers, and patients with complications of cancer chemotherapy and radiation treatment. LED's work by increasing the energy produced in the skins cells and accelerating tissue growth and regeneration.

When it comes to ageing of the skin LED treatment reverses signs of photo ageing, increases collagen production, regenerates tissues, reduces redness, breaks down pigmentation, heals acne and when combined with Impulse Micro-Currents lifts sagging tissue (think the dreaded waddle!).

How does the Soli-Tone 2500 compare to other LED treatments on the market?

LED's use specific wavelengths of Red, Blue, Green and Yellow light energy to stimulate the changes in your skin. Up to now, treatments (think Omnilux & MAX7) targeted individual problems (ie. They use one spectrum of light or "scan" the light over the skin). Research has shown that the best results are seen when all 4 light wavelengths are used in conjunction and the skin is exposed for a constant and set time, depending on the condition being treated.

Imagine if you tried to get a tan in a Solarium with just a few globes "scanning" back and forth across you - you would need 1000 times as many sessions to get your tan. The same applies when it comes to LED's that merely scan back and forth across the face and so the whole face is not constantly exposed to the LED and so light absorption is reduced.

The Soli-Tone 2500 is the only technology that uses high strength LED's (stronger than any other on the market), all wavelengths in a single treatment and allows for the light to be present for long enough AS WELL as using Micro-Currents to lift tissues. So it is the only technology that uses both light energy (LED's) and Micro-Currents to treat the epidermis, dermis and muscular regions of your skin.

How does the LumiFacial® treatment work?

The LumiFacial® targets common skin problems, such as blemishes, redness, and uneven pigmentation. During the treatment, a wrap-around facial panel is used. The Soli-Tone 2500 selects the specific Light energy required for your skin type and condition.

The treatment is completed in only 30 minutes. You can immediately return to your normal routine with no redness, peeling or necessity for time away from work. Treatment is very comfortable. Unlike Laser or IPL that rely on thermal damage, LumiLift® and LumiFacial® are non-invasive and work to heal - not damage. There is no risk of burns and is completely safe.

So, you get better results, quicker AND underlying muscles are lifted. There is no other technology available that uses LED's combined with Micro-Currents.

What exactly is the Lumilift® procedure?

The LumiLift® procedure is a revolutionary way of rejuvenating the skin and lifting sagging muscles. It is performed with two special electrodes that feature an advanced application of Micro-current and Light energy.

During the Lumilift treatment, a pulsating light is directed at the skins surface. Light passes through the layers of the skin and penetrates into the dermis. This powerful light stimulates the fibroblast cell to produce more collagen. The increased collagen smoothes and softens the appearance of wrinkles, reduces the size of pores, and improves the skin texture. At the same time, the electrodes emit Micro-currents to tone and lift sagging muscles. At the end of the treatment, an antioxidant product is infused into your skin using a special electrode to smooth away facial creases.

What are Micro-currents?

Micro-currents are non-invasive. They stimulate and resemble the biological currents that your body produces to tighten muscles. Extensive clinical research has shown that these currents are able to tone and lift the muscles. Combining Light Therapy with Micro-currents erases years from the appearance of your skin, without the use of drastic, costly and painful surgical procedures. Lumilift® will restore a natural and youthful appearance to your skin.

How many treatments do I need?

Age, lifestyle, and skin condition will influence the number of treatments needed. While clients usually report a noticeable improvement after just one treatment, a series is required to get the full benefits. Generally, depending on your skin and the severity of the conditions being treated, 6-10 treatments achieve optimal and long-lasting results.

What results can I expect?

The improvement can and will be quite dramatic. Most clients see immediate improvement. With a series of treatments, years are easily removed from the appearance of your skin. You will see younger looking skin, reduction in the amount of, and the depth of, fine lines and wrinkles, increased skin elasticity, and redefined and firmer facial contour. Pore size is reduced along with textural irregularities, acne, and pigmented blemishes. Have a look at the before and after photos on the front page.

This is achieved without any thermal damage to the skin (as with Lasers and IPL) and no down-time. Soli-Tone treatments are unique since they regenerate not only the face, but also the neck, décolletage and hands.

How long will the effects of the LumiLift® treatment last?

Improvements continue for at least 3 months after a single session. After completing a series, depending on the rate of the ageing process, changes can last over one year. A good skincare regime combined with a monthly maintenance treatment, is sufficient to enhance and maintain the results you obtained with the LumiLift® treatment initially.

How long does a LumiLift® or LumiFacial® treatment take?

The LumiFacial® takes only 20-30 minutes. When included with the LumiLift® you need to allow 40 - 50 minutes. You may also include a LumiLift® or LumiFacial® into your normal monthly maintenance facial as an add on treatment. As soon as your treatment is finished you may apply make-up and return to your normal everyday activities.